

HD 3406.2



HD 3406.2 **BENCH-TOP CONDUCTIVITY METER**

The HD3406.2 is a bench top instrument for electrochemical measures: conductivity and temperature.

The displayed data can be stored **(datalogger)** and can be transferred to PC or serial printer thanks to the multi-standard serial port RS232C and USB2.0 and software DeltaLog9 (Vers.2.0 and subsequent ones). The storing and printing parameters can be set from menu.

The HD3406.2 measures conductivity, liquid resistivity in liquids, total dissolved solids (TDS) and salinity using combined 4-ring and 2-ring conductivity/temperature probes. Temperature is measured by Pt100 or Pt1000 immersion, penetration or contact probes.

The probe calibration can be performed automatically in one or more of the $147\mu S$, $1413\mu S$, 12880μS or 111800μS/cm conductivity calibration solutions.

The display shows continually the temperature in °C or °F and one selectable parameter according to the connected probe type, i.e. in case of conductivity probe it is possible to select between χ or Ω or TDS or NaCl.

Other functions of this instrument include: Max, Min and Avg function, the Auto-HOLD function, the automatic turning off which can also be excluded.

The instruments have IP66 protection degree.



Technical characteristics HD3406.2 X, Ω, TDS, NaCl, °C/°F measurement

Instrument Dimensions (Length x Width x Height) Weight Materials Display

Operating conditions Working temperature Storage temperature Working relative humidity Protection degree

Power **Batteries** Autonomy (only batteries) Mains (cod. SWD10)

Security of memorized data

Selectable storage interval

Date and hour Accuracy

Serial interface RS232C Type Baud rate Data bit **Parity** Stop bit Flow Control Serial cable length

Selectable print interval

USB Interface Type

Common connections to all models Serial interface and USB Mains adapter (cod. SWD10)

Storage of measured values Type Quantity

Measurement connections Input conductivity Input for temperature probes complete with TP47 modules

Measurement of conductivity by instrument Measurement range (Kcell=0.01) Measurement range (Kcell=0.1) Measurement range (Kcell=1)

220x120x55mm 460g (complete with batteries) ABS, rubber 2x41/2 characters plus symbols visible area: 52x42mm

-5 ... 50°C -25 ... 65°C 0 ... 90% RH without condensation

3 batteries 1.5V type AA 100 hours with 1800mAh alkaline batteries Output mains adapter 100-240Vac/ 12Vdc-1A

Unlimited

1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

Schedule in real time 1min/month max drift

RS232C electrically isolated Can be set from 1200 to 38400 baud None

Xon/Xoff Max 15m

immediate or 1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour


1.1 - 2.0 electrically isolated

8-pole MiniDin connector 2-pole connector (positive at centre) 12Vdc/1A

2000 pages of 18 samples each 36,000 sets of measures made up of $[\chi - \Omega]$ or TDS or NaCl] and [°C- °F]

8-pole male DIN45326 connector 8-pole male DIN45326 connector

Resolution 0.000...1.999µS/cm $0.001 \mu S/cm$ $0.00...19.99\mu\text{S/cm}$ 0.01 uS/cm 0.0...199.9µS/cm 0.1μS/cm 1μS/cm 200...1999µS/cm 2.00...19.99mS/cm 0.01mS/cm 20.0...199.9mS/cm 0.1mS/cm



Range di misura (Kcell=10). 200...1999mS/cm 1mS/cm Accuracy (conductivity) ±0.5% ±1digit Measurement of resistivity by instrument Measurement range (Kcell=0.01) Up to $1G\Omega$ -cm Measurement range (Kcell=0.1) Up to $100M\Omega$ cm Measurement range (Kcell=1) 5.0...199.9Ω·cm 0.1Ω cm $200...999\Omega\text{-cm}$ $1\Omega\text{-cm}$ $1.00k...19.99k\Omega\text{-cm}$ $0.01k\Omega\cdot cm$ $20.0k...99.9k\Omega$ ·cm $0.1k\Omega\cdot cm$ 100k...999kΩ·cm 1kΩ.cm $1...10 M\Omega \cdot cm$ 1MΩ·cm Measurement range (Kcell=10) $0.5...5.0\Omega$ ·cm $0.1\Omega\cdot cm$ Accuracy (resistivity) $\pm 0.5\% \pm 1$ digit Measurement of total dissolved solids (with coefficient X/TDS=0.5) 0.005 ma/lMeasurement range (Kcell=0.01) 0.00...1.999mg/l Measurement range (Kcell=0.1) 0.00...19.99mg/l 0.05mg/l Measurement range (Kcell=1) 0.0...199.9 mg/l 0.5 mg/l 200...1999 mg/l 1 mg/l 2.00...19.99 g/ $0.01 \, g/l$ 20.0...99.9 g/l 0.1 g/l Measurement range (Kcell=10) 100...999 g/l 1 g/l Accuracy (total dissolved solids) ±0.5% ±1digit Measurement of salinity Measurement range / Resolution 0.000...1.999g/l 1mg/l 2.00...19.99g/l 10mg/l 20.0...199.9g/l 0.1g/IAccuracy (salinity) ±0.5% ±1diait Temperature measurement by instrument Measurement range Pt100 -50...+200°C Measurement range Pt1000 -50...+200°C Resolution 0.1°C ±0.25°C Accuracy Drift after 1 year 0.1°C/year

Automatic/manual temperature compensation

 $0...100^{\circ}\text{C with }\alpha_{\text{T}} = 0.00...4.00\%^{\circ}\text{C}$ Reference temperature $20^{\circ}\text{C or }25^{\circ}\text{C selectable from menu}$

Conversion factor χ/TDS 0.4...0.8

Cell constant K (cm⁻¹) 0.01 - 0.1 - 0.7 - 1.0 - 10.0

Standard solutions automatically detected (@25°C)

147μS/cm 1413μS/cm 12880μS/cm 111800μS/cm

(*) The resistivity measurement is obtained from the reciprocal of conductivity measurement. Close to the bottom of the scale, the indication of resistivity appears like reported in the table below:

K cell = 0.01 cm ⁻¹		K cell = 0.1 cm ⁻¹	
Conductivity (µS/cm)	Resistivity (MΩ·cm)	Conductivity (µS/cm)	Resistivity(M Ω ·cm)
0.001 μS/cm	1000 MΩ·cm	0.01 μS/cm	100 MΩ·cm
0.002 μS/cm	500 MΩ·cm	0.02 μS/cm	50 MΩ·cm
0.003 μS/cm	333 MΩ⋅cm	0.03 μS/cm	33 MΩ⋅cm
0.004 μS/cm	250 MΩ·cm	0.04 μS/cm	25 MΩ·cm
	•••		

ORDERING CODES

HD3406.2: The kit is composed of: instrument HD3406.2 datalogger, for measurement of conductivity - resistivity - TDS - salinity - temperature, 3 1.5V alkaline batteries, operating manual and DeltaLog9 version 2.0.

pH/mV electrodes, conductivity probes, dissolved oxygen probes, temperature probes, standard reference solutions for different measurement types, connection cables for pH electrodes with S7 connector, cables for data download to PC or printer have to be ordered separately.



χ







ACCESSORIES

HD2110CSNM: 8-pole connection cable Mini Din - Sub D 9-pole female for RS232C, for connection to PC without USB input.

HD2101/USB: Connection cable USB 2.0 connector type A - 8-pole Mini Din for connection to PC with USB input.

SWD10: Stabilized power supply at 100-240Vac/12Vdc-1A mains voltage.

HD40.1: Portable, serial input, 24 column thermal printer, 57mm paper width.

HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm.

HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.

TP47: Module for the connection of Pt100 4-wire and Pt1000 2-wire probes.

Combined conductivity and temperature probes

SP06T: Combined conductivity and temperature 4-electrode cell in Platinum, body in Pocan. Cell constant K = 0.7. Measurement range 5μ S/cm ...200mS/cm, 0...90°C.

SPT401.001: Combined conductivity and temperature 2- electrode cell in stainless steel AlSI 316. Cell constant K = 0.01. Measurement range 0.04 μ S/cm ...20 μ S/cm, 0...120°C. Measurement in closed-cell..

SPT01G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 0.1. Measurement range $0.1\mu S/cm...500\mu S/cm, 0...80^{\circ}C$.

SPT1 $\vec{\mathbf{g}}$: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 1. Measurement range 10 μ S/cm ...10mS/cm, 0...80°C.

SPT10G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 10. Measurement range 500μS/cm ...200mS/cm, 0...80°C.

Electrode characteristics at page 402

Standard conductivity calibration solutions

HD8747: Standard calibration solution 0.001mol/l equal to 147μ S/cm @25°C - 200cc. HD8714: Standard calibration solution 0.01mol/l equal to 1413μ S/cm @25°C - 200cc. HD8712: Standard calibration solution 0.1mol/l equal to 12880μ S/cm @25°C - 200cc. HD87111: Standard calibration solution 1mol/l equal to 111800μ S/cm @25°C - 200cc.

Temperature probes complete with TP47 module

TP47.100: Direct 4 wires Pt100 sensor immersion probe. Probe's stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47.1000: Pt1000 sensor immersion probe. Probe's stem Ø 3mm, length 230mm. Connection cable 2 wires with connector, length 2 m.

TP87.100: Pt100 sensor immersion probe. Probe's stem Ø 3mm, length 70mm. 4 wire connection cable with connector, length 1 m.

TP87.1000: Pt1000 sensor immersion probe. Probe's stem Ø 3mm, length 70mm. 2 wire connection cable with connector, length 1 m.

Accessorie

TP47: Module for the connection of Pt100 4-wire and Pt1000 2-wire probes.

mg/l